

# OPERATING INSTRUCTIONS

## Center of Mass Demonstrator No. 74836

### 1. Introduction

The Center of Mass Demonstrator is designed to illustrate the position of the centers of mass of differently shaped objects.

### 2. Description

The Center of Mass Demonstrator consists of four black geometric shapes: an equilateral triangle, a sector of a circle, a half-circle, and an irregular shape. Also included are an adapter from which the shapes are suspended, a plumb line and plumb, and several pieces of chalk.

### 3. Experiment

Determine the center of gravity or mass of each geometric shape by suspending the shapes from a ringstand. For example, the equilateral triangle is hung successively from the predrilled holes in each apex. Each time, the chalk line is attached at the point of suspension and snapped against the figure, leaving a vertical line of chalk. The center of mass is where the lines intersect. This can then be checked by punching a hole where the lines intersect, suspending the shape from the hole, and noticing that it remains in equilibrium independent of orientation about the point of suspension. In the case of the irregular shape, it is found that the center of mass lies outside the body since the chalk lines do not intersect on the shape.

### 4. Discussion

The above experiment proves that the center of gravity of an equilateral triangle lies at the intersection of the medians. The center of gravity of a sector of a circle is at  $\bar{x} = (2R \sin \theta) / 3\theta$  where  $\theta$  is the semivertex angle, in this case  $30^\circ$ . Accordingly, for this sample  $\bar{x} = 2R / \pi$  (see Figure 1).




Figure 1

The semicircle is a sector of total angle  $\pi$ . Accordingly,  $\bar{x}$  becomes  $4R / 3\pi$ . The irregular sample yields this interesting observation: *The center of gravity of a body does not have to lie where there is matter.* So, the plumb lines intersect where there is nothing — the center of gravity of a doughnut, for example is in the middle of the hole.



## **6. Maintenance**

The Center of Gravity Demonstrator needs no special maintenance. Should any difficulty develop with this apparatus, please contact Central Scientific Company giving details of the problem. To ensure better service, please do not return any item to Central Scientific until we have sent written authorization.

## **7. Copyright Notice**

The Center of Gravity Demonstrator operating instructions are copyrighted and all rights reserved. Permission is granted to all non-profit educational institutions to make as many copies of these instructions as they like as long as it is for the sole purpose of teaching students. Reproduction by anyone for any other reason is prohibited.

Revised 2/92